



Production and Manipulation of Radioactive Ion Beams

Almost all activities to be discussed can naturally be adapted for EURISOL

Iain Moore

University of Jyväskylä, Finland

Rare isotope beam preparation and manipulation

- Beam purity *mass sep., chemical selectivity*
- Emittance *cooling*
- Energy spread *cooling & trapping*
- Temporal structure *pulsing/bunching*
- Ionic/atomic state *optical pumping*
- Charge state *ionization/breeding*


A personal viewpoint from recent developments!

RIBs: from high to low energy

Gas catchers: universal devices capable of extracting recoils from essentially all production mechanisms into a low-energy beam.

GOAL

High-precision experiments with exotic nuclei almost at rest
(production by projectile fragmentation / fission)




- universal (all but He) and fast production (~5-50 ms)
- high selectivity (inherent via in-flight separation)
- high efficiency (though element dependent)
- could be used in connection with EURISOL !

Novel high intensity gas catchers


CARIBU



- low-energy and re-accelerated RIBs (masses, decay and laser spectroscopy)
- 80 mCi ^{252}Cf source, 70 isotopes for PAC expts.
- **340 mCi source (September 2012)**
- Yield with new source 4-5 \times previous source:
~20000 ^{141}Cs /s at 27+ out of the charge breeder,
~10000 ^{141}Cs /s at 4.5 MeV/u at Gammasphere



^{252}Cf




G. Savard et al., *Int. J. of Mass. Spec.* 251 (2006) 252
http://www.phy.anl.gov/atlas/caribu/ATLAS_Cf_upgrade.pdf

UHV cells to cryogenic catchers ($\sim 70\text{K}$)





Advantages

- Ultra-pure helium
 - Ideal for ion survival
 - No formation of molecules/adducts
- 2+ charge state
 - shorter extraction times



A cryogenic gas cell for SHIPTRAP



Transport efficiency of α -decay recoil ions in a closed chamber

P. Dendooven et al., NIMA 558 (2006) 580

A. Saastamoinen, I.M. et al., NIMA 685 (2012) 70

- Fully assembled, rf ready
- Cooled to 30K, at 50 mbar room temp equivalent
- ^{223}Ra recoil source tests
- Ready for on-line use

S. Eliseev, C Droese et al.

S. Eliseev et al., Nucl. Instr. and Meth. B 266 (2008) 4475–4477





Cryogenic stopping cell design for S-FRS



JUSTUS-LIEBIG-
UNIVERSITÄT
GIESSEN


Kernfysisch Versneller Instituut
university of
 groningen





M. Ranjan et al., Europhys. Lett. 96 (2011) 52001

M.P. Reiter, Master Thesis, Justus-Liebig-Universität Gießen (2011)

On-line experiments: October 2011, July 2012



Stopping cell performance





Preliminary


	^{223}Th
Stopping efficiency	$(27 \pm 3)\%$
Survival and extraction efficiency	$(43 \pm 9)\%$
Total efficiency	$(12 \pm 2)\%$

- Almost $2\times$ > density than other rf- cells
- Mean extraction time from centre of cell $\sim 25 \text{ ms}$
- Intensity limitation not expected for most Super-FRS beams

Selective resonant laser ionization



Gas cell selectivity and spectroscopy



Separation of stopping and laser ionization volumes improves:

- laser ionization efficiency for high cyclotron beam current
- Ion selectivity

Production of ^{94}Rh
Selectivity:

[Laser(on)/Laser(off)]


Ion Collector OFF = 450

Ion Collector ON = 2200

- In-gas-cell spectroscopy of neutron deficient Cu isotopes (*T. Cocolios et al., PRL 103 (2009) 102501*)
- Investigate possibility of performing similar measurements near $N=Z$ line around $A=100$ (LISOL and JYFL)

Heavy Element Laser Ionization Spectroscopy

An alternative approach towards laser spectroscopy of heavy elements



- Expected production rates:

$^{208}\text{Pb}(^{48}\text{Ca}, 2n)^{254}\text{No}$: ~ 150 pps (10 μA primary beam intensity)



- Strategy

In-gas cell laser ionization spectroscopy (broadband – 5 GHz):

In-gas jet laser ionization spectroscopy (narrow band – 200 MHz)


Hot cavity laser ion source trap

“To achieve the highest selectivity in RIB production”




F. Schwellnus et al., Rev. Sci. Instrum. 81 (2010) 02A515.

LIST on-line at ISOLDE (2012)



- LIST was successfully tested with UC_x target
 - > No loss of performance over 5 days
- Suppression of Na-, Al-, K-, Fr-, U-isotopes studied
 - > Suppression factors varied from 100 to 1000
- Laser ionization of radioactive Mg and Po in LIST mode


Repeller



Transducer box

rf connector


LIST inside target





Courtesy of Daniel Fink

LIST: suppression of Fr; ionization of Po

In-source laser spectroscopy: a Po survey



²¹⁶Po, first attempt in 2007



Courtesy of Thomas Cocolios & Daniel Fink


Low energy ion beam coolers

Rare isotopes often produced with non-ideal properties



- Duty cycle not matched to experiment (eg: ion traps)
- Energy spread too high (eg: laser spectroscopy)
- Transverse emittance too high (eg: RIB facilities)
- “debuncher” for CW formation, P. Delahaye next talk !

More than a dozen RFQ coolers in operation/construction

- Efficiency $\sim 40\text{-}80\%$ with beam deceleration
- Energy spread $< 1\text{eV}$
- Transverse emittance $\sim \text{few } \pi \text{ mm mrad}$
- Bunch widths $\leq \text{few } \mu\text{s}$




TITAN's RFQ buncher (TRIUMF)




E. Mané et al., PRC 84 (2011) 024303

Optical manipulation of atomic states




- $J=0 \rightarrow J=1$ gives μ , Q , $\delta\langle r^2 \rangle$ but not I
- Access to more accessible/efficient transitions
- New elements to study (eg. transition metals)



F.C. Charwood et al., Phys. Lett. B 674 (2009) 23
B. Cheal et al., Phys. Rev. Lett. 102 (2009) 222501

Future extension: polarization of beams




High intensity beam cooler - SHIRAC

Spiral II High Intensity Radiofrequency Cooler (SHIRAC)

- Designed to manipulate the high intensity RIBs from SPIRALII target to HRS

Requirements:


- Current up to $1 \mu\text{A}$
- 60% transmission ($A=132$)
- Emittance 1π mm mrad
- Energy spread $\sim 2\text{eV}$
- No bunching, DC mode only



Special RF system designed for the high intensity


- 6.5 MHz at 8kV
- Tunable and harmonic

5 loops





SHIRAC - current status

Surface ionization source: Cs⁺ at 5 keV



- Emittance $\sim 1.5 \pi$ mm mrad for $1 \mu\text{A}$
- Nuclear environment adaptation underway


Work of PhD student Ramzi Boussaid



MR-TOF Mass Spectrometry at ISOLTRAP



- Multi-reflection time-of-flight mass separator at the ISOLTRAP experiment



Support Penning-trap mass spectrometry

Mass measurements get access to beams with minute production rates, tens of ms half-lives and high suppression of contaminants due to:

- auxiliary isobaric purification
- fast purification
- accumulating of ion of interest

MR-TOF plus detector as stand-alone system


The MR-TOF offers a way to separate background for direct single-ion detection using MCP (time scale: tens of ms)

- direct mass measurements
- in-source spectroscopy


R. N. Wolf et al., Nucl. Instr. and Meth. A 686, 82-90 (2012)




Recent results at ISOLTRAP



In-source spectroscopy for hyperfine structure



Direct mass measurements of $^{52-54}\text{Ca}$




See poster of S. Kreim (ISOLTRAP)

S. Kreim et al., INTC-P-299, IS 518 (2011) INTC-P-317, IS 532 (2011)

Time of flight / arb. units



MR-TOF-MS and the ion catcher at GSI



Mass Resolving Power

600,000

Mass Measurement Accuracy

$\sim 10^{-7}$

Measurement Duration

~ 10 ms


Transmission efficiency

up to 70%

- Broadband mode for characterization of the cryogenic ion catcher (diagnostics tool)
- High resolution mode for mass measurement (^{213}Rn with half-life 19.5 ms)


Collinear Resonance Ionization Spectroscopy

CRIS: “Combining the high resolution nature of collinear beams with the high sensitivity of in-source spectroscopy”



Decay spectroscopy station:
(Laser-assisted decay measurements)

Courtesy of Kieran Flanagan




Key points:

- Efficiency achieved 1:60
- Background suppression $\sim 1:500\,000$ (pressure limited, expect 1:1e6)
- $202-206, 218, 219, 229, 231\text{Fr}$

Laser-assisted decay spectroscopy

- Using the HFS, decay spectroscopy can be made on pure isomeric (or ground state) beams
- Assisted α -, β -, and γ -decay spectroscopy



- Important advancements continually made at all facilities to be able to handle the ever-increasing improvements in primary/secondary beam intensities:
 - *purity, energy spread, temporal structure, charge state, ionic/atomic state....*
- Gas catchers becoming more common and finding new applications:
 - *versatile, universal*
- Cooler-bunchers are standard laboratory equipment:
 - *High efficiency, universal, easy to build*
 - Next frontier approaches – *high intensity* devices
- Commissioning of MR-TOF-MS devices complements Penning trap techniques and can be used in a variety of modes (stand-alone devices, diagnostic tools etc)
- Novel in-source/in-jet spectroscopy demonstrated at both ISOL and gas cell-type facilities; plans for future facilities.
- ❖ Not discussed: application of storage rings, Penning traps and so on!!!
- ❖ Many of these techniques can be applied at a future EURISOL facility 😊