

Università di Pisa  
LM Materials and Nanotechnology - a.a. 2016/17

## Spectroscopy of Nanomaterials II sem – part 10

Version 0, May 2017  
Francesco Fuso, francesco.fuso@unipi.it  
<http://www.df.unipi.it/~fuso/dida>

**Optical near-fields and related techniques:  
seeing close, seeing small**

# OUTLOOK

We well know how diffraction can affect optical spectroscopy of nanomaterials, ruling the ultimate diffraction achievable in conventional (general purpose) configurations

Diffraction is strongly involved in near-fields, which can be considered as created by strongly diffracting apertures

Quite surprising is to discover that near-fields can be used to obtain sub-diffraction spatial resolution in virtually all (almost all) spectroscopy configurations

Furthermore, being near-fields relevant close to small-sized apertures, emitters, scatterers, they play a rather ubiquitous role in nano-optics and nanophotonics

Today's menu:


- Appetizers: near-field cooked à la Fourier, à la Heisenberg, à la Maxwell, and other variants
- Main course of tapered fibers, used in scanning probe microscopy, garnished with several examples of optical spectroscopy of different materials
- Dessert: tip enhanced spectroscopies, Raman flavoured

# NEAR-FIELD I

## Basic idea of the near-field:


To stay **near** the emitter, where “near” means at a (small) fraction of the wavelength

*Typically, for ordinary emitters near-field is indistinguishable from conventional (propagating) radiation...*



*...but...*

If the emitter is, e.g., a small, sub-wavelength sized, aperture, near-field “prevails”







# ALREADY SEEN

The concept of near-field and the related consequences are ubiquitous in the nanophotonics world

*For instance:*

*Field enhancement in localized plasmon resonances occurs in the near-field*

*Refraction in a LHM produces “amplification” in the near-field (actually non discussed)*



## NEAR-FIELD II

According to Fourier optics:

$$g_{out}(x, y, z) = \int_{-\infty}^{+\infty} \int_{-\infty}^{+\infty} T(k_x, k_y) G_{in}(k_x, k_y, z) e^{i2\pi(k_x x + k_y y)} dk_x dk_y$$

$$G_{in}(k_x, k_y, z) = \int_{-\infty}^{+\infty} \int_{-\infty}^{+\infty} g_{in}(x, y) e^{-i2\pi(k_x x + k_y y)} e^{i2\pi k_z z} dx dy$$

Due to the small aperture size, the transfer function  $T$  mimics the transfer function of a **low-pass filter** with very high cut-off frequency (spatial frequency):

$$T(k_x, k_y) \simeq 1 \quad \text{for} \quad -u < k_{x,y} < u$$

with  $u \sim 1/a$

**Reconstruction of sub-wavelength details in the image is possible according to Fourier**

On the other hand

$$G_{in}(k_x, k_y, z) = \int_{-a/2}^{+a/2} \int_{-a/2}^{+a/2} g_{in}(x, y) e^{-i2\pi(k_x x + k_y y)} e^{i2\pi k_z z} dx dy \approx f(k_x, k_y) e^{-i2\pi(k_x a + k_y a)} e^{i2\pi k_z z}$$

**The near-field gets an evanescent, non-propagating character**

$$k = \frac{2\pi}{\lambda} = \sqrt{k_x^2 + k_y^2 + k_z^2}$$

$$\rightarrow k_z = \sqrt{\left(\frac{2\pi}{\lambda}\right)^2 - (k_x^2 + k_y^2)}$$

for  $k_x, k_y \rightarrow u$ ,  $k_z \in \text{Im}$

## NEAR-FIELD III

### Hole (or object) size $a > \lambda$

- low spatial frequencies
- root  $> 0$
- $k_z$  real
- progressive wave along  $z$

$$g_{out}(x, y, z) \sim f(k_x, k_y) e^{i(k_x x + k_y y)} e^{ik_z z}$$


### Hole (or object) size $a < \lambda$

- high spatial frequencies
- root  $< 0$
- $k_z$  imaginary
- evanescent wave along  $z$

$$g_{out}(x, y, z) \sim f(k_x, k_y) e^{-\alpha z}$$

### Note:

in the far-field there is still a propagating wave, strongly diffracted and attenuated because of the small aperture size (*Fraunhofer diffraction*)



# HEISENBERG, DIFFRACTION AND NEAR-FIELD

Diffraction can be seen as a sort of consequence of Heisenberg's (uncertainty) principle: this is the content of a conceptual experiment (Heisenberg's microscope)

**Propagating waves:**

$k_x$  are real and  $|k_x| \leq k = 2\pi/\lambda$



Heisenberg's principle:  $k_x \geq 2\pi/\Delta x$

$$\Delta x \geq \lambda$$

(actual parameters give the Abbe's limit)

*... but ...*

**In non-propagating (e.g., evanescent) waves:**


$k_x$  can be imaginary, and, e.g.:  $|k_x| \geq k = 2\pi/\lambda$



***The Heisenberg's principle is no longer ruling the ultimate resolution!***


Extreme diffraction means wavevectors for the field leaving the aperture (or a nanosized object) gets a large in-plane component

→ ***Propagation is somehow frustrated***



**In some sense, near-field enables using ultra-diffractive fields for beating the diffraction limit**

# CONVERTING AN EVANESCENT FIELD






The evanescent field at a surface can be “converted” into an ordinary (propagating) e.m. wave through a sort of **photon tunneling** process across a vacuum gap

- *near-field at a surface can be “captured” by a suitable probe*
- *near-field can couple to a surface and gets a propagating character*

Or, if you have a fluorescent object, **the fluorescence can be excited by the near-field** (it is the time-dependence, i.e., the frequency that must be “resonant”) and collected in the far-field

Or, if you have a scatterer (a nanosized object) **near-field can be scattered in the far-field**

# NEAR-FIELD AND MAXWELL EQS



$$\mathbf{E} = \frac{1}{4\pi\epsilon_0} k^2 (\mathbf{n} \times \mathbf{p}) \times \mathbf{n} \frac{e^{ikr}}{r} + [3\mathbf{n}(\mathbf{n} \cdot \mathbf{p}) - \mathbf{p}] \left( \frac{1}{r^3} - \frac{ik}{r^2} \right) e^{ikr}$$

$$\mathbf{H} = \frac{ck^2}{4\pi} (\mathbf{n} \times \mathbf{p}) \frac{e^{ikr}}{r} \left( 1 - \frac{1}{ikr} \right)$$

In the near-field ( $r \rightarrow 0$ ,  $kr \ll 1$ )

**“Static” dipole behavior**  
**→ no propagation**

$$\mathbf{E} = \frac{1}{4\pi\epsilon_0} [3\mathbf{n}(\mathbf{n} \cdot \mathbf{p}) - \mathbf{p}] \frac{1}{r^3}$$

$$\mathbf{H} = \frac{i\omega}{4\pi} (\mathbf{n} \times \mathbf{p}) \frac{1}{r^2}$$

# NUMERICAL CALCULATIONS

The actual field distribution is perturbed by the presence of the imaged object  
 ↓  
 Numerical simulations (e.g., multiple multipole, FDTD, etc) needed




Fig. 3 – Intensity maps for the electric field ( $E^2$  and components  $E_x^2$ ,  $E_y^2$ , and  $E_z^2$ ) at a distance of  $z = 0.3a$  from the tip aperture of radius  $a$  created by (a) a surface charge distribution as illustrated in fig. 2, and (b) a planar Bethe aperture (Bethe-Bouwkamp solution). All images show an area of  $4a$  by  $4a$ , and the images in each row share the same color bar. The two images on the left correspond to those that would be obtained using a pointlike, scalar detector for the electric field (e.g. an idealized, infinitely small fluorescent nanosphere). The images of the last three columns correspond to fluorescence images predicted for pointlike vector detectors for the electric field (e.g. single, fluorescent molecules) oriented along the  $x$ ,  $y$ , and  $z$  direction, respectively.

**Actual space distribution of the near-field can differ from the ideal case because of interaction with objects and of border effects**

EUROPHYSICS LETTERS

1 April 2004

*Europhys. Lett.*, **66** (1), pp. 41–47 (2004)

DOI: 10.1209/epl/i2003-10138-7

The optical near-field of an aperture tip

A. DREZET, M. J. NASSE, S. HUANT and J. C. WOHL(\*)



# LOCALIZATION OF THE E.M. FIELD

In practical terms near-field enables a spatial localization of the e.m. field in a size **ruled by the aperture**

*Typically, localization is on a  $\lambda/5$ - $\lambda/10$  scale*

- *Ideal resolving power  $\sim \lambda/5$ - $\lambda/10$*
- *Resolving power no longer ruled by diffraction!!*

The Scanning Near-field Optical Microscope (SNOM, aka NSOM) is based on such concept



# BRIEF HISTORY

1928/1932

E.H. Syngé proposes the idea of using a small aperture to image a surface with sub-wavelength resolution using optical light. For the small opening, he suggests using either a pinhole in a metal plate or a quartz cone that is coated with a metal except for at the tip. He discusses his theories with A. Einstein, who helps him develop his ideas. [E.H. Syngé, "A suggested method for extending the microscopic resolution into the ultramicroscopic region" *Phil. Mag.* 6, 356 (1928); E.H. Syngé, "An application of piezoelectricity to microscopy", *Phil. Mag.*, 13, 297 (1932)].

1956


J.A. O'Keefe, a mathematician, proposes the concept of Near-Field Microscopy without knowing about Syngé's earlier papers. However, he recognizes the practical difficulties of near field microscopy and writes the following about his proposal: "The realization of this proposal is rather remote, because of the difficulty providing for relative motion between the pinhole and the object, when the object must be brought so close to the pinhole." [J.A. O'Keefe, "Resolving power of visible light", *J. of the Opt. Soc. of America*, 46, 350 (1956)]. In the same year, Baez performs an experiment that acoustically demonstrates the principle of near field imaging. At a frequency of 2.4 kHz ( $\bar{\lambda} = 14$  cm), he shows that an object (his finger) smaller than the wavelength of the sound can be resolved.

1972

E.A. Ash and G. Nichols demonstrate  $\bar{\lambda} / 60$  resolution in a scanning near field microwave microscope using 3 cm radiation. [E.A. Ash and G. Nichols, "Super-resolution aperture scanning microscope", *Nature* 237, 510 (1972)].

1984

The first papers on the application of NSOM appear. These papers are the first to show that NSOM is a practical possibility, spurring the growth of this new scientific field. [A. Lewis, M. Isaacson, A. Harootunian and A. Murray, *Ultramicroscopy* 13, 227 (1984); D.W. Pohl, W. Denk and M. Lanz, *APL* 44, 651 (1984)].



Sketch of the first proposal contained in a letter from Syngé to Einstein (1928)

## The Nobel Prize in Chemistry 2014



Photo: A. Mahmoud  
Eric Betzig  
Prize share: 1/3




Photo: A. Mahmoud  
Stefan W. Hell  
Prize share: 1/3



Photo: A. Mahmoud  
William E. Moerner  
Prize share: 1/3

Eric Betzig was among the discoverers of SNOM (not very satisfied, though)

# ENHANCED RESOLUTION



SNOM and confocal PL images of butterfly-shaped monolayered MoS<sub>2</sub> prepared on an ITO-coated glass substrate






**Spatial resolution is actually enhanced compared to confocal microscopy**

***In actual operating conditions, the resolving power can reach to few tens of nm range***

Lee, et al., Nanoscale (2015)

# SNOM POWER (VIRTUALLY)



Virtually, all configurations of conventional microscopy can be translated into a near-field version, e.g., **collection**, **illumination**, scattering, fluorescence, optical activity, etc.





**But, in any case, surface properties are probed**

# SNOM PROBES (APERTURE)

Hollow pyramids (cantilevered) or tapered optical fibers are used as SNOM probes




**We will restrict here to tapered optical fibers probes**




Typical aperture diameter  
50-100 nm

# TAPERED OPTICAL FIBERS

## Conventional probe fabrication (*heat and pull*)




## Chemical etching



The final result is a conical (tapered) optical fiber, with the taper metal-coated in order to prevent losses, and an apical (uncoated) hole playing the role of the aperture

# TAPER BEHAVIOR

A whole slew of scanning near field optical microscopes (SNOMs) have been developed by researchers over the past 15 years. These microscopes smash the diffraction limit of far field microscopes, potentially achieving resolution an order of magnitude better than a standard confocal microscope [Hecht 2000]. One such SNOM is the scanning tunneling optical microscope. This microscope uses a sharp glass tip to locally frustrate total internal reflection below a surface, indirectly imaging features on that surface at high spatial resolution.




Aperture based SNOMs are more common and more practical. [Hecht 2000] They produce an evanescent field by forcing light through a small aperture (see figure). The evanescent field locally illuminates the sample. Once free of the aperture, the field is no longer evanescent, and it expands in the far field to be picked up by a detector. To achieve high resolution, the aperture must be small, and close to the sample surface so the field is tightly confined when it interacts with the sample. Full analysis of the field-sample interaction of a SNOM is a difficult or impossible undertaking, but the data produced can yield important and detailed information about a sample surface.

Suggested reading: Hecht et al., J. Chem Phys. 112, 7761 (2000).



Note: metal layer (typ Cr, Ni) can absorb radiation → power entering the fiber cannot exceed the mW range!

Note: probe “throughput” (i.e., ratio between output/input power) is quite low for fiber probes,  $\sim 10^{-5} - 10^{-8}$ , but near field **intensity** can be large enough for sensitive analyses

A simplified ray optics picture in the tapered region.  
Light can both:  
- be back reflected  
- (partially) absorbed by the metal layer




# LIMITATIONS (A FEW)



L. Novotny and D. W. Pohl, in *Photons and Local Probes*, NATO ASI Series E, p.21-33, Kluwer Academic, 1995.

In a real, tapered fiber probe, understanding the actual field distribution inside and right outside the probe is difficult, since:

- It is not an aperture in a thin opaque screen
- The taper angle and all (polarization-dependent) processes occurring in the taper affect the field
- Details of the aperture shape can be relevant



L. Novotny and D. W. Pohl, in *Photons and Local Probes*, NATO ASI Series E, p.21-33, Kluwer Academic, 1995.

# TIP-TO-SURFACE DISTANCE

A distinctive feature of near-field is that its amplitude rapidly drops to zero within a range  $\sim a$



During the scan, the probe tip must be kept “close” to the surface (typ at a distance  $< 10$  nm)

“Constant gap” operation is strictly required for the SNOM images to be reliable




A method to continuously monitor tip/sample distance is needed

A feedback acting on the vertical piezo displacement is used (as, e.g., in any SPM)



***A topography image is simultaneously acquired during each scan, with a lateral resolution depending on the probe size, typ in the 100 nm range***





If the tip is kept in **longitudinal** oscillation, the oscillation amplitude depends on the distance due to shear-forces (mostly associated with friction of the air layers between tip and surface)

## Notes:



Oscillation amplitude must be small (typ  $\sim 1$  nm) to prevent resolution loss

Oscillation at resonance frequency is required to get maximum sensitivity

# SHEAR-FORCES AND TOPOGRAPHY




Fiber tip is glued onto a *tuning fork* (a quartz diapason acting as a mechanical oscillator and showing capabilities to measure the oscillation amplitude)



- A dithering piezoelectric transducer keeps the probe tip in oscillation along a direction **parallel** to the surface
- Oscillation amplitude is monitored by a *tuning fork*
- When the distance gets smaller (typ., below 10 nm), the oscillation is **damped** (and phase is changed) due to **shear-forces** involving many effects (e.g., viscous interaction of the air layer between tip and sample)
- Similar to AFM in tapping mode, but for the oscillation direction, the relevant distance and the involved mechanisms

**A topography map is acquired simulatneously and in the same region of the optical maps**

# ACCESSIBLE CONTRAST MECHANISMS



Virtually, there are no limitations to the range of optical features that can be measured at the local scale, i.e., to the optical spectroscopies that can be realized in the near-field

Note, however, that only surface phenomena can be analyzed (unless complicated models for data interpretation are deployed)


A range of possible configurations is feasible in order to operate different nanoscale spectroscopies

A selection of examples is given here in the following

# EXAMPLES IN PLASMONICS

## Sample:

pattern of lithographed gold film on glass, with a grating connected to a waveguide/concentrator featuring a tapered end




Far-field excitation and near-field collection are accomplished at different positions  
 → propagation can be investigated at the local scale

## Goal:

to investigate the propagation of **surface plasmon** modes along plasmonic waveguides/concentrators


## Method:

**collection mode SNOM**



# SPP PROPAGATION

SNOM map




- Surface plasmon modes collected in the near-field down to the tapered end
- Concentrator effect at the tapered end is evident


✓ Confirmation of plasmonic waveguide behavior

✓ Complete characterization accomplished

SNOM map



- Interference between propagating and counterpropagating plasmon modes well seen in flat-end waveguides




[J. Berthelot, et al., JOSA B (2012)]

# LSP HOT-SPOTS

## Sample:

gold nanowire arrays made by ion beam sputtering at grazing incidence (*produced at Univ. Genova by F. Buatier de Mongeot*)


Topography map




## Goals:

- to determine the spatial distribution of e.m. field at the local scale
- to identify the **hot-spots** due to the collective plasmonic resonances in the inhomogeneous gold nanoparticle arrangement

Near-field Red





Near-field Green



Polarization orthogonal to NW axis

**Hot-spots clearly detected and spatially mapped**  
*(only with proper wavelength and polarization of the far-field excitation)*



Polarization parallel to NW axis

[C. D'Andrea, et al., J. Phys. Chem. C (2014);  
M. D'Acunto, et al., Beilstein J (2017)]

# SMART POLYMERS I

## Samples:




host-guest systems consisting of dye molecules dispersed in a polymer matrix  
(produced by A. Pucci, G. Ruggeri at DCCI)

## Goal:

to identify emission domains at the local scale and to relate them with the matrix conditions


## Method:

Local excitation of the sample by the near-field, photoluminescence collected in the far-field



The degree of molecular aggregation, depending on concentration, matrix structure and/on the application of mechanical stress, affects the photoluminescence properties


# SMART POLYMERS II



[G. Iasilli, et al., Macromol. Chem. Phys. (2014)]


# SMART POLYMERS III

## R-pery in Low Lin Dens PE (LLDPE)




## Thermochromic effect (at the macro-scale)

Under UV illumination  
R-Pery 0.1% wt. in LLDPE



## Macroscopic emission spectrum




## Goal:


to identify emission domains at the local scale and to investigate their properties as a function of temperature

## Method:

Local excitation of the sample in the near-field, photoluminescence collected in the far-field with simultaneous two color detection




# SMART POLYMERS IV



- ✓ Red emitting domains more intense in the “valleys”
- ✓ Green emitting domains more intense on the “hills”
- ✓ Excitation at high near-field intensity leads to local spectral modifications  
(not shown here)

# OPTICAL ACTIVITY



Optical activity, i.e., the response of a material to **polarized radiation**, is a key point in many applications and fundamental issues

Its **analysis and measurement** at the local scale is a challenging task


✓ **Basic idea:**

*to excite the sample in the near-field with a **polarization modulated radiation***

*to apply demodulation techniques in order to extract the physical quantities of interest*

The method has been applied to a variety of properties such as, birefringence, linear and circular dichroism


# POLARIZATION MODULATION



Core of the setup:  
photoelastic modulator (Hinds Instrument PEM-100)

$$\vec{E}_{PEM} = \frac{E_0}{\sqrt{2}} (\hat{e}_x + \exp(i\varphi(t))\hat{e}_y) \quad \varphi(t) = A\cos(2\pi ft)$$

Produced polarization states (at different times)





At the SNOM probe input, laser polarization is periodically manipulated (circular-elliptical-linear)

Circular polarization achieved twice in a modulation period  $T$  (left- and right-handed)


Assuming as input a linearly polarized radiation at 45 degrees,  $A = \pi/2$ ,  $T = 2\pi/f$  with  $f = 50$  kHz

# EXAMPLE I : ELECTROSPUN FIBERS



Sample:  
MEH-PPV polymer prepared in  
micro- and nano-fibers through  
electro-spinning



Example of topography maps



# MOTIVATIONS




## Model:

- inhomogeneous structure across the fiber width due to spontaneous stretching processes
- fiber core made of strongly aligned polymer molecules


***expected linear dichroism in the core***

**Polarization analysis with sufficient spatial resolution needed to assess the model**

# RESULTS



Line profiles of topography, absorption, dichroic ratio




**The model is substantially confirmed**

**Engineering the optical properties of the fiber is feasible**



[A.Camposo, et al., Nano Lett. (2013);  
A.Camposo, et al., Macromol. (2014)]

# CHIRALITY AND CIRCULAR DICHROISM I




Molecular chirality:  
describes a molecular structure whose mirror image cannot be superimposed


carbon  
nanotubes



amino acids



enzymes




DNA

**Chirality is extremely frequent in biomolecular building blocks (homochirality in nature is an intriguing issue)**



# CHIRALITY AND CIRCULAR DICHROISM II

Chirality is inherently associated with circular dichroism (CD)

**CD spectroscopy is a powerful tool to investigate chirality in the liquid phase**




For the sample of our interest (TPPS<sub>3</sub> aggregates) in solution:



**CD very hard to be observed in the solid state**

## EXAMPLE II : PORPHYRIN NANOSTRUCTURES

- ✓ Chirality can clearly stem from mutual arrangement of different achiral molecules
- ✓ In the case of the TPPS<sub>3</sub> porphyrin system, a self-assembly molecular aggregation process in the liquid phase can occur




**Hypothesis: chirality originates from molecules self-assembled in nano-helices**


**Need for high resolution optical methods able to probe individual, deposited, nanostructures**

# CIRCULAR DICHROISM AT THE SOLID STATE


Circular dichroic ratio



Linear dichroic ratio



Demodulated SNOM optical maps



Linear dichroism related to the elongated (anisotropic) structure


Circular dichroism intrinsic to the material nanostructure

[F.Tantussi, et al., Nanoscale (2014);  
G.Lazzini, et al., SPIE Proc., (2016)]

# APERTURELESS AND TIP-ENHANCED

A variant of SNOM is based on aperture-less probes

Such probes are available, **semiconducting** or, more frequently, metal-made, with nanometer apical size (the same technology of STM or AFM probes)



A metal tip (similar to STM probes) is illuminated by a far-field

→ Scattering produces a near-field localized at the tip


→ Assuming prevalent interaction with the (higher intensity) near-field, spatial resolution is obtained

*(unfortunately, discrimination between near-field excited and far-field excited features is not always simple)*


# APERTURELESS RESOLUTION

Apertureless-SNOM is mostly used with fluorescent samples, thanks to the straightforward (bandpass filter based) discrimination of illumination and fluorescence signals

The extremely small size of the tip ensures excellent spatial resolution




Spatial distribution of electric field close to the tip



SNOM image of a single dye molecule anchored on a substrate



# TERS EXAMPLES





A typical example of TERS is with Carbon Nanotubes, showing particularly intense Raman bands related to their structure

Note that typical size of single walls CNT is on the nanometer range

Sheremet. et al., Chemnitz (2014)

# TERS+STM



Molecular islands of H<sub>2</sub>TBPP deposited on Ag

# CONCLUSIONS

- ✓ Near-field entails the e.m. field in proximity to small structures, and therefore plays a crucial role in the optical interaction between closely spaced, nanosized objects
- ✓ Near-field can be produced by a sort of extreme diffraction from nanosized apertures and its properties in imaging can be understood by looking at the problem from different standpoints:
  - (i) Fourier optics
  - (ii) Heisenberg principle
  - (iii) Maxwell equations
- ✓ A (virtually) extremely powerful spectroscopy of nanomaterials can be built based on near-fields
- ✓ However, technical and practical limitations exist to run such spectroscopies, and results, while being useful, are often controversial and non completely satisfying

*In any case, the meaning of making optical spectroscopy of nanomaterials should be finally clarified!!!*

## FURTHER READING

For a comprehensive reference on near-field and nano-optics

A. Zayats and D. Richards, Nano-optics and near-field optical microscopy, Artech House, Boston (2009).

A paper on SNOM from the NT-MDT company:

<http://www.azonano.com/article.aspx?ArticleID=2250>

A comprehensive list of references on SNOM can be found at

<http://zeiss-campus.magnet.fsu.edu/referencelibrary/superresolution/nsom.html>

A review on SNOM, authored by Novotny and Stranick, can be found at

<http://www.optics.rochester.edu/workgroups/novotny/muri03/publications/stranick06a.pdf>

A (rather old) review on SNOM applications to materials, authored by Dunn, is

<http://pubs.acs.org/doi/abs/10.1021/cr980130e>